

# EcoVac-200

## Portable Suction Pump



The Merridian EcoVac-200 Portable Suction Pump is designed for low to medium suction applications, for use in nursing home, homecare and clinic environments. It combines economy with efficiency.

### FEATURES

- Powered by durable oil-less diaphragm pump with strong suction pressure.
- Equipped with automatic overflow device to prevent overflowing when bottle is full with liquid.
- Equipped with effective bacterial filter.
- Compact, durable plastic collection bottle, and handle on housing for portability.

### SPECIFICATIONS

**Max. Vacuum**

560 mmHg

**Sound Level**

≤65dB(A)

**Power Supply**

AC 220 ± 10%, 50Hz ± 2%

**Net Weight**

3.9 Kg

**Displacement**

≥18L/min

**Storage Bottle**

Durable Plastic 1000ml

**Dimension (L x W x H)**

28cm x 19.5cm x 22cm

**Certificates**

ISO 13495, CE

## EcoVac-200DC

(AC/DC Version)



Power Supply:  
AC & 12V DC  
from vehicle cigarette  
or built-in Battery Pack.

Merridian  
Health Care