

EcoVac-200

Portable Suction Pump



The Merridian EcoVac-200 Portable Suction Pump is designed for low to medium suction applications, for use in nursing home, homecare and clinic environments. It combines economy with efficiency.

FEATURES

- Powered by durable oil-less diaphragm pump with strong suction pressure.
- Equipped with automatic overflow device to prevent overflowing when bottle is full with liquid.
- Equipped with effective bacterial filter.
- Compact, durable plastic collection bottle, and handle on housing for portability.

SPECIFICATIONS

Max. Vacuum
560 mmHg

Sound Level
≤65dB(A)

Power Supply
AC 220 ± 10%, 50Hz ± 2%

Net Weight
3.9 Kg

Displacement
≥18L/min

Storage Bottle
Durable Plastic 1000ml

Dimension (L x W x H)
28cm x 19.5cm x 22cm

Certificates
ISO 13495, CE

EcoVac-200DC

(AC/DC Version)



Power Supply:
AC & 12V DC
from vehicle cigarette
or built-in Battery Pack.