

# JPD-500C Pulse Oximeter



With MINI-USB Port



## Features

- ❖ Integrated with SPO<sub>2</sub> probe and processing display module.
- ❖ Small in volume, light in weight and convenient in carrying.
- ❖ Operation of the product is simple, low power consumption.
- ❖ SPO<sub>2</sub> and PR display.
- ❖ Pulse rate value display, bar graph display.
- ❖ Low-voltage indication: low-voltage indicator appears before working abnormally which is due to low-voltage.
- ❖ Automatically power off function: when the device is under the state of measuring interface.
- ❖ Various color of cover can be selected.
- ❖ Auto shut off :16s
- ❖ OLED display, Adjustable brightness
- ❖ SPO<sub>2</sub> : 35-100%
- ❖ Pulse Ratio : 25-250BPM
- ❖ Resolution: 1% for SPO<sub>2</sub>, 1BPM for pulse ratio
- ❖ Accuracy: ±1% (70%-100%), unspecified (<70%) for SPO<sub>2</sub>, ±1BPM or ±1% (select larger) for pulse ratio
- ❖ Power: 1.5V (AAA size) alkaline batteries × 2
- ❖ Size: 58.0×31.6×32.9mm

## Packaging information

Unit dimension:

G/W(KG):0.09

N/W(KG):0.06

Size:83×77×50mm

Carton dimension: (100PCS/carton)

G/W(KG):10

N/W(KG):8.5

Size:420×390×280mm