

Mannual Resuscitator



Fortune[®] Manual Resuscitator Features:

- Equipped with pressure relief valve at 40 cmH₂O with overridden mechanism.
- 3 different materials of Silicone, Rubber and PVC for reuse or disposable selections.
- 3 sizes are available: infant, child and adult.
- ISO standard connectors.

Model	Infant	Child	Adult	Description
Silicone Resuscitator	SR-Infant	SR-Child	SR-Adult	Reusable
Rubber Resuscitator	RR-Infant	RR-Child	RR-Adult	Disposable
PVC Resuscitator	PR-Infant	PR-Child	PR-Adult	Disposable

Silicone Resuscitator:

			Silicone mask #0/ #1/ #2 Patient valve w/ pressure relief valve Silicone bag 240 ml Reservoir valve Oxygen reservoir 600 ml
			Silicone mask #3 Patient valve w/ pressure relief valve Silicone bag 500 ml Reservoir valve Oxygen reservoir 2600 ml
			Silicone mask #5 Patient valve w/ pressure relief valve Silicone bag 1600 ml Reservoir valve Oxygen reservoir 2600 ml

