

Valueklave 1730MK

The Speed and Quality you Demand
with the Value you Deserve



tuttnauer

Sterilization & Infection Control Products

The Valueklave1730MK will satisfy all of your basic sterilization needs at an economical price. High speed and reliability are built into every Valueklave1730MK



Valueklave1730MK Specifications

| | |
|-------------------------------------|---|
| Chamber Dimensions Ø x L In (mm) | 6.7"x13.4" (170 x 340) |
| Volume | 7.5L/2US gal |
| Overall Dimensions DxWxH In (mm) | 17.9"x17.4"x12" (455 x 440 x 305) |
| Standard Cassettes | 2 half |
| Tray Dimensions DxWxH In (mm) | 4.7"x11.6"x 0.8" (120 x 295 x 20) |
| No. of Trays | 3 |
| Standard Unwrapped Cycle Time | Cold: 12 mins Hot: 9 mins |
| Power Requirement | 120V 1Ph 50/60Hz 11.3A 230V 1Ph 50/60Hz 5.9A |
| Power Consumption | 1350W |
| Shipping Weight lbs (kg) | 55 (25) |

Features:

- Fast Cycle Time. 9 minutes from a hot start (Standard Unwrapped Cycle time).
- Automatic shut off at the end of both the sterilization and dry cycles
- A long life electro polished chamber and door.
- Double safety locking device prevents door from opening while chamber is pressurized.
- Drain valve is located on the front, allowing for quick and easy draining of water reservoir.
- Dual safety thermostat to protect against overheating.
- International certifications.



www.tuttnauer.com

tuttnauer
Sterilization & Infection Control Products

Tuttnauer Europe b.v.
Paardeweide 36, 4824 EH
P.O.B. 7191 4800 GD
Breda, The Netherlands
Tel: +31 (0) 765423510
Fax: +31 (0) 765423540
E-mail: info@tuttnauer.nl

Tuttnauer USA Co.
25 Power Drive
Hauppauge, NY, 11788, USA
Tel: (800) 624 5836
(631) 737 4850
Fax: (631) 737 0720
E-mail: info@tuttnauer.com