



*




german|engineering

BF 105
Body Complete

beurer
wellbeing

BF 105 Body Complete Diagnostic bathroom scale



Complete body analysis: Accurate upper and lower body measurement with 8 electrodes



Innovative networking between PC, smartphone and scale
With free "beurer HealthManager" app + software for downloading
Transfer using Bluetooth® low energy technology or USB



Detailed body analysis with graphic evaluation



Determination of body weight, body fat, body water, muscle percentage and bone mass

Separate indication of upper and lower body measurement for body fat and muscle percentage possible

Basic turnover and active turnover in kcal

Removable XXL display (suitable for wall mounting) with blue LCD illumination

Time and date

Intuitive user guide in 5 languages (D/GB/F/ES/I)

Additionally with USB cable and "beurer HealthManager" software

download for free or app to download

With brushed stainless steel electrodes

Change over to kg/lb/st

Capacity: 180 kg

Graduation:

- 0.1% for body fat, body water and muscle ratio

- 100 g for body weight and bone mass

5 activity levels

10 user memories with 30 measurement spaces each + guest mode

Overload indication

Button-On-Technology

Automatic switch-off

Incl. 6 x 1.5V batteries AA

Product measurements (LxWxH):

Scale: 33 x 33 x 4 cm Display: 21.2 x 7.7 x 3.8 cm

Product weight:

Scale: approx. 1.73 kg Display: approx. 260 g

5 year guarantee

Sales unit: 2 / Shipping carton: -

EAN no.: 4211125 74812 8

Item no.: 748.12

* tested online application "beurer HealthManager" Cloud

0617 Subject to errors and changes



*




german|engineering

BF 105
Body Complete

beurer
wellbeing



CONNECTION



Innovative networking between PC, smartphone and scale - keep track of your measurements and body data at home and on the move. Easily synchronised with your smartphone via Bluetooth®



Syncs with

- iOS ≥ 8.0
- Android™ ≥ 4.4
- Bluetooth® 4.0



List of supported devices

0617 Subject to errors and changes