



HK 72 Akku Mobiler Wärmegürtel

Bis zu 4 Stunden mobile Wärme ohne Steckdose - für maximale Bewegungsfreiheit

Leistungsstarker Lithium-Ionen-Akku mit LED-Anzeigen für Ladefunktion und Temperaturstufen

Mit Ladegerät, Akku Ladezeit ca. 2-3 Stunden

Körpergerechte Form für Bauch und Rücken

Aus dünnem Material, kann unauffällig unter der Kleidung getragen werden

Anschmiegsam und hautsympathisch

Breites Gummiband mit Klettverschluss - angenehmer Tragekomfort

4 beleuchtete Temperaturstufen

Abschaltautomatik nach ca. 90 Minuten

Schnellheizung

Maschinenwaschbar bis 30°C

4 Watt

Netzteilspannung:
Input: 100-240V, 50-60 Hz, 0,6 A
Output: 8,4V DC, 1,2 A

Produktmaße (LxB): 115 x 18 cm

3 Jahre Garantie

VE: 3 / Transportkarton: -

EAN-Nr.: 4211125 22591 9

Art.-Nr.: 225.91

HK 72 Akku Mobile Heat Belt

Up to 4 hours mobile heat independent from the mains - for maximum mobility

Powerful lithium-ion-battery pack with LED displays for charge function and temperature settings

With charger, charging time approx. 2-3 hours

Economic shape for abdomen and back

Thin material, comfortable to wear under clothing

Cosy and skin friendly

Wide rubber strap with hook and loop fastener - comfortable to wear

4 illuminated temperature settings

Automatic switch-off after approx. 90 minutes

Rapid heat up

Mashine washable up to 30°C

4 watt

Power supply voltage:
Input: 100-240V, 50-60 Hz, 0.6 A
Output: 8.4V DC, 1.2 A

Product measurements (LxW): 115 x 18 cm

3 year guarantee

Sales unit: 3 / Shipping carton: -

EAN no.: 4211125 22591 9

Item no.: 225.91