




HM200P Portable Height Stadiometer

The classic portable height stadiometer can be disassembled into smaller pieces for easy transportation using the included carry case, perfect for measurement on the go!





FEATURES

- Wall stabilizers for use against wall for maximum stability
- Convenient and easy to setup and dismantle, with transport using carry box
- Large floor plate ensures stability
- Markings printed on each side for easy viewing
- Light weight and detachable parts make installation and measurement simple



- Portable carry box

SPECIFICATIONS

| | |
|---------------------------|---|
| Model | HM200P Portable Height Stadiometer |
| Measuring Range | 14 – 205 cm / 5.5 – 81 inch |
| Graduation | 1 mm, 1/16 inch |
| Units | mm, inch |
| Standard Accessory | Support plate x 1 / Head stopper x 1 / Base platform x 1 / Height rod x 4 |
| Weight | 2.0 kg / 4.4 lb (Carry Case 0.6 kg / 1.3 lb) |
| Dimensions | 350(W) x 400(D) x 2430(H) mm Carry Case: 645(W) x 365(D) x 125(H) mm |