

PULSE OXIMETER

MD300

C26



PULSE OXIMETER

Product Features:

- Display SpO2 waveform (fill or line)
- OLED display
- Six kinds of display mode
- Switch off automatically
- Brightness adjustable
- Low voltage indicator
- Low power consumption
- 2 AAA alkaline or rechargeable batteries



超思
CHOICEMMED

Technical Specifications

Patient Range

Adult, Pediatric and Neonatal

SpO₂

Range 0%-100%
Resolution 1%
Accuracy 80%-100%: $\pm 2\%$
70%-80%: $\pm 3\%$
0%-69%: unspecified

Heart (Pulse) Rate

Range 0-255 bpm
Resolution 1bpm
Accuracy ± 2 bpm or $\pm 2\%$

Alarm

Alarm Low power indicator

Display

Type OLED
Parameters SpO₂, PR, PI, SpO₂ waveform
Mode 6 display modes
Brightness adjustable 10 levels

Environmental

Operating Temperature 5-40°C
Storage Temperature -20°C - 55°C
Operating Humidity 15%RH - 80%RH
Storage Humidity 10%RH - 80%RH


Mechanical

Size 58*32*34mm
Weight: 50g


Battery

Type 2AAA alkaline batteries
Operation time about 30 hours for normal operation


Overview the data in any direction



Mode1



Mode2



Mode3



Mode4



Mode5



Mode6



AUTHORIZATION AGENT