MS3500

Infant Scale



Functions

•Hold •Unit switch (kg/lb)

•Tare •Removable Tray

Model	MS3500 Infant Scale
Capacity	20 kg
Graduation	5g < 10 kg > 10g
Overall Dimensions	560(D) x 350(W) x 120(H) mm
Tray Dimensions	560(D) x 290(D) x 65(H) mm
Platform Dimensions	325(D) x 310(D) x 32(H) mm
Device Weight	2.6 kg
Power Supply	9V battery or power adapter
Accessories	Height Rod (HM80D, HM80M), Carrying Bag

The MS3500 digital infant scale features sloped edges to safely and comfortably hold infants during weight measurement. Highly accurate, the scale can track even slight changes in infants' weight.

The scale can be used to weigh toddlers as well after removing the measurement tray, and features a hold function to give accurate readings even for moving subjects.



Measure toddlers with measurement tray removed



HM80D Digital Infant Height Rod



HM80M Mechanical Infant Height Rod