



# MS 4910

- Capacity: 300 kg
- Graduation: 100 g
- Functions: BMI, Print, Weight & Height measurement, Hold, Tare
- Accessories: HM201M / HM200D Height Stadiometer  
TP2100 / TP2110 Thermal Printer




- Bubble level indicator
- Adjustable rubber feet to adjust level
- Anti slip rubber mat
- Print readout of Weight, Height and BMI via thermal printer
- Compatible with digital and manual height stadiometers



Displays Weight, Height, and BMI on dual LCD dscreens



## • Features



## + Accessories

