

MS6111

Double-sided Medical Floor Scale

The MS6111 is an affordable, accurate, professional medical scale, designed with a low profile platform to make it easier to step on and off. The double-sided LCD screen makes it possible for results to be viewed easily by either the patient or the operator. When used with the HM202P Portable Stadiometer, weight and height measurement can be conducted simultaneously, saving valuable time and effort.



Double-sided indicator for easy viewing



Compatible with HM202P Stadiometer – no cumbersome base attachment needed

Functions

- BMI
- Bubble Level
- Print
- Adjustable Feet
- Hold
- Anti-slip Mat
- Tare
- Unit

Specifications

Capacity	250 kg
Graduation	10 g < 60 kg > 100g
Overall Dimensions	345(W) x 350(D) x 45(H) mm
Device Weight	2.5 kg
Accessories	AR-2781 Carry bag, HM202P Height Stadiometer TP2110 Thermal Printer



AR-2781 Carry Bag