

SL 60 Snore mask


**german
engineering**


Termination of snoring via
vibration alarm



Adjustable pulse intensity with
preview



Precise snore analyses by
"beurer SleepQuiet" app via *Bluetooth*[®]



Gentle snore therapy to reduce snoring for a good night sleep

Innovative real-time detection of snoring as it extracts the user's snoring voice and bone conduction vibration signal

Immediate termination of snoring via vibration alarm at the forehead. Due to body reaction the throat muscles become firm and it may change the lying position.

Comfortable sleeping mask
- Machine washable up to 30 °C
- Adjustable length

With free app „beurer SleepQuiet“ for:
- Precise snoring analysis
- Display and graphical representation of data (day/week/year)
- Determination of snoring volume
- Snore diary
- Adjustable pulse intensity with preview

Compatible with iOS 9.0 and Android™ 4.4 with *Bluetooth*[®] 4.0

Rechargeable battery
- Charging time: 2 hours
- Runtime: 2 - 3 nights in case of heavy snoring

Incl. micro USB charging cable

Product measurements (LxWxH): Controller: 4.3 x 1.9 x 1.2 cm
Mask: 23.5 x 9.0 x 3.4 cm


Product weight: approx. 48 g

3 year guarantee

Sales unit: 4 / Shipping carton: 4 x 4

EAN no.: 4211125 67769 5

Item no.: 677.69



The innovative snore mask SL 60 acts gently on the sleeper. If it discovers a snoring noise by tone detection and structure-borne sound, it sends vibration pulses to interrupt the snoring. Due to body reaction the throat muscles become firm and the airways open. Breathing is calmed and the snoring decreases. It may also change your lying position. Have a peaceful night!

Homescreen for quick overview and navigation

Colourful snoring indicator:



Snoring sounds per hour

0 - 19
20 - 49
50 - 99
≥ 100

Colour
●
●
●
●

Snore level
normal
light
middle
heavy




Diary function:
Documentation of influencing factors as alcohol, stress, etc.



Graphical overview:
Day, week, year



Settings:
Individualisation of vibration alarm (with preview), sensitivity, etc.



Bluetooth®
Transfer using
Bluetooth® low
energy technology



List of supported smartphones



System requirements:
iOS ≥ 9.0
Android™ ≥ 4.4
Bluetooth® 4.0

*The app you download may differ from the illustrations.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners. Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Google Play and the Google Play logo are trademarks of Google LLC. Android is a trademark of Google LLC.